

NOVA FORM (NV)

MD

WINN CARE FRANCE, Miniac
ZONE ACTIPOLE
2 RUE DE SAINT-COULBAN
35540 MINIAC-MORVAN - FRANCE

WINN CARE FRANCE, Nîmes
200 rue Charles Tellier
Actiparc de Grézan
30034 NIMES Cedex 1 - FRANCE

IFU_NOVA

IFU_NOVA

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1. INSTRUCTIONS

Purpose of the device

Prevention of pressure ulcers. Handicap compensation.

Indications

The visco-elastic foam mattresses are bed supports which help to prevent pressure ulcers in patients who are at moderate to high risk for pressure ulcers (who have been evaluated using a validated scale and clinical judgment), who are up standing during the day, and bedridden for more than 15 hours in every 24 hour period.

- › One or more stage 1 and/or 2 pressure ulcers, not located on the weight bearing areas, or with the possibility of support exclusion (using a localized pressure relief system), patient standing or not during the day.
- › One stage 3 or 4 pressure ulcer not located on the weight bearing areas, with the possibility of support exclusion (using a localized pressure relief system),
- › One stage 1 or 2 pressure ulcer in the weight bearing areas, patient standing or not during the day, with a positioning system and intervention by medical assistants 3 times per day for verification of the installation and for turning.

Contra-indications

- › Periodic repositioning is not possible,
- › Patients at very high risk of pressure ulcers,
- › Patients with stage 3 or 4 pressure ulcers in the weight bearing areas,
- › Patients whose weight is higher than that which the mattress can support.

Patient and user target groups

Adult patients who are hospitalized, institutionalized, or at home, presenting with one or more pressure ulcers and/or who are at risk of developing pressure sores due to any temporary or permanent alterations in their health. These individuals are cared for by health professionals, with the help of assistants if required.


Indicate the undesirable side-effects

All serious incidents which are related to the device must be notified to the manufacturer and the competent authority of the member state in which the user and/or patient resides.

Inform the competent authority if you believe, or if you have reason to believe, that the device presents a serious risk, or that it has been tampered with.

2. COMPOSITION OF THE MEDICAL DEVICE

- 1 Polyurethane foam mattress.
- 2 Cover with zipper (non-removable for the PROMUST HIP).
- 3 Identification label on the outside of the cover.
- 4 Instructions for use (paper or electronic).

3. CLINICAL BENEFITS, PERFORMANCE, MECHANISM OF ACTION

Performance characteristics of the device

The visco-elastic foam mattresses allow for pressure relief in supporting areas to improve blood circulation and thus reduce the risk of pressure ulcers in the weight bearing areas of the body, through immersion and envelopment.

Expected clinical benefits

Maintenance of tissue oxygenation in the anatomical areas which are in contact with the supporting surface, by reducing the pressure applied to the cutaneous and sub-cutaneous tissue.

Information for health professionals

Observe the condition of the patient's skin which is in contact with the support structure several times a day. Verify that the mattress is correctly installed.

4. USAGE INSTRUCTIONS

User training

User training must be provided by people who are trained and validated by the relevant economic operators. The choice of support must take into account the risk of pressure ulcers, as per the Braden scale, or any other validated scale, the individual's weight, and its use by the caregivers.

Installation of the device

The mattress is delivered with its zipped cover. The product is ready to be installed. The visco-elastic foam surface must be placed in contact with the patient's body. The cover must be dry. Place the mattress on the bed in the proper direction, with the cover's flaps facing down. Lie the patient down.

MATTRESS DELIVERED COMPRESSED

Foot/head reversible mattress

⌚ 2 mn → 10 mn *

* This time may be increased if the storage temperature is lower than 15°C and/or in line with the storage period (0 to 6 months).

The return to a normal visual appearance of the mattress under normal usage temperature conditions may take up to 24 hours, but this does not present any risk for his use. Small creases may be visible on the lower surface of the CIC cover after unrolling. They present no risk to the patient and do not affect the efficiency of the final product.

MATTRESS DELIVERED FLAT

Foot/head reversible mattress

Preventive maintenance

- › Annually carry out a visual inspection of the state of the foam: visible subsidence of the material, and a non-equal and very slow return to the foam's original state are signs of ageing which compromise its pressure ulcers prevention properties.
- › Once a year, check the external surface of the cover, by using a light source from inside in order to check for holes and/or tears.

Cleaning and disinfection

- › Between each patient.
- › Wash the cover in a machine at 90°C.
- › Do not iron the cover.
- › Can be used with a bleach solution of up to 5,000 ppm.

Information regarding suitable procedures for re-use

- › Each time that it is re-used, the product must be physically and bacteriologically cleaned by performing the procedures given in the technical manual (which can be downloaded at www.winnicare.com).
- › Replace the cover if the appearance of its surface changes.

5. WARNINGS AND PRECAUTIONS

Required measures

A support by itself is not enough to prevent pressure ulcers. Other preventative measures are also required.

- › Change positions at least once every 4 hours.
- › Maintain skin hygiene and avoid maceration.
- › In case of incontinence, regularly change any protective items.
- › Once a day, observe the condition of the skin, or ensure that it is observed.
- › Ensure that the patient's food intake is sufficient and appropriate.
- › Regularly drink enough fluids.
- › A minimal distance of 22 cm between the top of the bed rail and the top of the mattress with the patient on it is advised.
- › Limit large distances between the body and the support by as much as possible, with the exception of the bed support sheet, clothing, or adult diapers.
- › Ensure that there are no foreign bodies such as: tubing, crumbs, grease, etc.
- › Reinforced surveillance in subjects at risk of falling (patient profile: invalid, agitated, not lucid).

Circumstances under which a user should consult with a health professional

If one of the required measures cannot be followed, you must inform your doctor or nurse as soon as possible. Inform your doctor or nurse as soon as possible of any abnormal events such as fever, pain, reddening or whitening of the weight-bearing areas (head, shoulder, back, hip, shoulder blade, pelvis, hip etc.)

Precautions

- › Bone and/or muscle trauma which has not been stabilized and which is in contact with the support.
- › First days following surgery for a pressure ulcer (skin graft or flap).
- › Patient being followed-up at home, without the possibility of receiving help from medical assistants.
- › Also verify the condition of the exposed part of the skin each time that care is provided, and during positional changes.
- › COPD-type bronchopulmonary disorders, with severe reduction in the tone of the respiratory muscles.
- › Joint stiffening and tendon retractions.
- › Fractures of the pelvis or spine.

Warnings

- › Not compatible with external cardiac massage, unless a board is placed between the thorax and the upper surface of the mattress.
- › Re-evaluate the patient's risk of entrapment in the non-moving parts of the medical bed associated with "therapeutic mattress and accessories and articulated positions of the bed frame" in accordance with the EN 60601-2-52 standard for adults.
- › The use of dry hot water bottles heated in the microwave in direct contact with the surface of the mattress is prohibited.

Lifespan information

The device can be reused as long as its initial performance is maintained throughout its life cycle. Its estimated lifespan is 6 years.

Accessories

Only the mattress protection accessory supplied by the manufacturer ensures the safety and performance of the complete product.

6. STORAGE/HANDLING/ DISPOSAL

Specific conditions for storage and/or handling

- › Store away from light, heat, and humidity.

Disposal of the product

Do not throw the product out into non-dedicated locations in nature. Use the recycling facilities available in your country.

Pressure (max) in mmHg before ageing, using the maximum validated weight

Max. user weight	Max. user pressure		
	Head	Buttocks	Heel
160 kg	95	77	85

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