
Axtair One®Plus

Axtair Automorpho®Plus

Axtair Axensor® AT12/AT15/AT20

Axtair XXL®

## Technical manual

CE

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

## ALARMS AND REPAIR

## Using the “info” key (Only on the Axensor line)

 A continuous press on this button when the compressor is running provides information about the last alarm recorded and the need for servicing.

Alarm codes: Description		Meaning
	Lock indicator steady	System problem Support
	Lock indicator flashing	Compressor system problem
	Dynamic indicator flashing	No pressure
	Treatment indicator steady	Initial inflation impossible
	Treatment indicator flashing	Inflation impossible
	Static indicator flashing	Over-pressure
	Dynamic and static indicator flashing	Pneumatic problem

The following table shows all the alarms for the device, their triggering conditions, their priority as well as possible repair.

Description	Trigger		Possible repair
	Condition	Time frame	
Alarm indicator flashing and mains indicator off  + sound signal	Electrical fault: the compressor is not powered by the mains	0 min	1 - Re-connect the compressor 2 - Check the status of the electricity network 3 - Replace the fuse
Alarm indicator flashing and mains indicator on  + sound signal	Compressor system problem	0 to 1 min	1 - Check that the mattress is properly connected to the compressor 2 - Check that the CPR valve is well sealed 3 - Check the leak-tightness of the mattress
	The compressor measures a nil pressure when pumping	1 min	
	The compressor measures an over-pressure	1 min	
	Inflation impossible	10 min	
	The compressor detects a pneumatic problem	1 min	
Alarm indicator steady  + sound signal	Initial inflation impossible	40 min	1 - Check that the mattress is properly connected to the compressor 2 - Check the status of the AXENSOR® box 3 - Check that the CPR valve is well sealed 4 - Check the leak-tightness of the mattress
	Support system problem	0 min	

For any other repair solution, contact your retailer or if you are trained, use the AirCare solution

**MAINTENANCE/DISINFECTION****COMPRESSORS:**

- **INDICATIONS:** Use a cloth soaked in a detergent product solution or CE\* labelled surface detergent/disinfectant at the use concentrations recommended by the manufacturer. Adhere to the remanence times.

- **CONTRAINDICATIONS:** 

**It is essential to disconnect the power cord connected to the compressor from the mains power outlet. The "mains" indicator must be off.**

- Do not use a high-pressure jet to clean an AXTAIR line compressor
- Do not spray liquid to avoid any damage to the compressor by penetration of liquid inside the box.

- **PRECAUTIONS:** 

Do not use abrasive cleaning products such as industrial detergents, acetone type solvents, ether and colouring products (iodized alcohol, potassium permanganate, silver nitrate etc.). Do not allow abrasive materials such as steel wool or "scotch brite" pads.

**MATTRESS :**

Depending on the type of mattress, there are several cleaning options available:

	Manual Disinfection 	Domestic washing machine 	Machine Barrier Aseptic 	REHAWASH Process 
Type of disinfection	Manual Disinfection	Covers Only	Complete mattress	Complete mattress
Mattress ONEP*	✓	✓	✓	✓
Mattress AUTOP*	✓	✓	✓	✓
Mattress AXENSOR	✓	✓		✓
Mattress XXL	✓	✓		

\* marketed from 01/11/22

All mattresses in the Axtair range are designed with an pneumatic assembly (Cells + Bundle + Maintaining Sheet) that can be easily removed from the base for easy cleaning and machine washing of the textile components.

For **AXENSOR® mattresses, the AXENSOR® pneumatic unit and control box** must be manually disinfected: use a cloth soaked in a detergent solution or CE\* labelled surface detergent/disinfectant at the use concentrations recommended by the manufacturer (Adhere to the remanence times) or follow the washing protocol of the REHAWASH process.

The **covers and the bases** can be machine washed in the following conditions:

	Wash in water, maximum temperature 90°C, reduced mechanical action, rinsing at decreasing temperature, reduced spinning.
	Bleaching possible, chlorination at 5000 ppm allowed.
	Ironing excluded.
	Dry cleaning not allowed, use of solvent-based stain remover not allowed.
	Drum drying authorized, moderate temperatures

**Aseptic Barrier Machine Washing:** ONEP and AUTOP mattresses marketed since 01/11/2022 are compatible with barrier machine disinfection.

For more information, please refer to the protocol sheet F039, presented in the following pages.

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FICHE 039

23/02/2023

## AXTAIR MATTRESSES – WASHING WITH ASEPTIC BARRIER MACHINES

## AIM

To wash the Axtair One Plus and Automorpho Plus mattresses thoroughly in an aseptic barrier machine. Effective washing ensures microbiological cleanliness of the mattress, thus preventing transmission of germs from one patient to another.

## EPI

## WHEN

Between each patient or before the product is stored.


## EQUIPMENT

- ▶ Professional aseptic barrier washing machine with disinfectant detergent,
- ▶ Professional tumble dryer,
- ▶ Washing net.

## PROTOCOL/PROCEDURE

axtair one<sup>®</sup>  
plus

1. Open the mattress cover completely and remove the air cell positioning sheet, held in place by two zips.
2. Remove the foam base, located in the compartment underneath the cells, from the mattress.
3. Close the air compressor connector to the mattress with the transport cap.
4. Loosely roll up the cell positioning sheet with the cells facing inwards. Place this roll inside the washing net and close with the tie provided.
5. Place the washing net and the mattress cover in the washing machine.
6. After the wash and spin cycle is complete place the washing net and mattress cover into the tumble dryer.
7. After complete drying reassemble the mattress and check that it is working correctly by inflating it with the rapid pump. Ensure all air cells are correctly aligned during this operation.
8. Repack the mattress in a clean transit bag for storage and delivery.

axtair automorpho<sup>®</sup>  
plus

1. Unzip the top cover.
2. Remove the air cell positioning sheet held in place by two zips.
3. Remove the foam base, located in the compartment between the anti slip base and the air cells, from the mattress.
4. Close the air compressor connector to the mattress with the transport cap.
5. Ensure that the heel relief valves are in the open position (*discharge not activated*).
6. Loosely roll up the cell positioning sheet with the cells facing inwards. Place this roll inside the washing net and close with the tie provided.
7. Place the washing net and the mattress cover in the washing machine.
8. After the wash and spin cycle is complete place the washing net and mattress cover into the tumble dryer.
9. After complete drying reassemble the mattress and check that it is working correctly by inflating it with the rapid pump. Ensure all air cells are correctly aligned during this operation. cette opération du bon positionnement des cellules.
10. Repack the mattress in a clean transit bag for storage and delivery.

## SETTINGS<sup>1</sup>

### Washing machine settings

- ▶ 60°C < wash temperature < 65°C (*provides thermal disinfection without damage to the mattress components*)
- ▶ Spin speed < 750 RPM
- ▶ Time : At 65°C the wash time must be at least 25minutes. Minimum wash time can be calculated using the following formula :

$$\sum m_i (T_i - 55) > 250$$

$\sum m_i$  = number of minutes during which the temperature in °C is maintained

$T_i$  = wash temperature in °C

The total number of mattresses that can be washed together depends on the capacity of the machine. As an example four mattresses can be washed in a 22kg capacity machine.

### Dryer settings

- ▶ 50°C < wash temperature < 65°C
- ▶ Time : Minimum drying time at 65°C estimated to be 48 minutes.


### IMPORTANT

- ▶ The tumble drying is an essential part of the process to ensure good microbial cleanliness of the mattress (*all parts of the mattress must be completely dry before reassembly*).
- ▶ The Air cell positioning sheet must remain in the washing net during the tumble drying process.
- ▶ Loading of the machine must comply with manufactures recommendations in terms of efficiency and cycle times.

<sup>1</sup> Maximum spin speed tested. It may be necessary to reduce spin speed depending on the machine used. These tests were carried out using a FAGOR machine.

## MAINTENANCE DIAGRAM AND RECOMMENDED DISINFECTION

\* Neutral detergent (household washing product)

\*\* Wide spectrum cleaning disinfectant product, compliant with standards Bactericide EN1040, EN13727+A2 ; Fungicide EN1275, EN13624 ; Sporicide EN14347 ; Virucide EN 14476+A2

## PREVENTIVE MAINTENANCE

### SERVICING FREQUENCY:

It is recommended to service the compressors and mattresses at least every 2 years of use for the following product lines:

- AXTAIR ONE PLUS®
- AXTAIR AUTOMORPHO AUTO PLUS®
- AXTAIR AUTOMORPHO AXENSOR®
- AXTAIR XXL®

### MATTRESS:

Top cover: visually inspect for damage and stains on the stitch side (white). If damage or stains are found, the cover must be replaced.

AXTAIR line mattresses do not contain any wear parts. A simple check of leak-tightness is recommended at the time of servicing.

Steps for checking mattress leak-tightness:

1. Inflate the mattress to a pressure between 85mBar and 100mBar.
2. Wait for the pressure to stabilize for approximately 30 seconds
3. Record the exact pressure of the mattress (P0)
4. Wait 1 minute and record the new pressure (P1)
5. The mattress is deemed to be compliant with the leak-tightness test if the difference between P0 and P1 (P0-P1) is less than or equal to 5 mBar.

In the event of non-conformity, it is necessary to identify the failing components (cells, connection, CPR valve etc.) and to replace them with components provided by the manufacturer. After replacement a further leak-tightness test should be performed.

### COMPRESSOR:

The filter of the compressor must be replaced once per year or more frequently depending on environmental conditions (dust, smoke, etc.).

This filter is on the back of the compressor behind a transparent card which shows its level of fouling. (The original filters are white)

AXTAIR compressors contain wearing parts that need to be replaced and a sensor that needs to be checked for calibration.

### SOLUTION AIRCARE

The Aircare solution was developed for the preventive and curative maintenance of these systems. Access to the solution requires:

- Following a specific one-day training session. Training module: WS/FORM/05. The detailed programme and dates for this training are available on [winnicare.fr](http://winnicare.fr) on the training tab. This training is also available as e-learning.
- Downloading specific software running in a Windows environment.
- The purchase of a connection kit containing the control tools.

The Aircare solution authorizes preventive and curative maintenance operations on the compressors and mattresses. The software is used for:

- quick troubleshooting of the checked system,
- the validation of servicing operations with removal of the maintenance key,
- repair validation,
- the traceability of all completed operations.

It is possible to find out if Axensor compressor servicing is required by a continuous press on the  button when the compressor is running.

Coding for servicing requirement: Description		Meaning
	Hourglass flashing quickly	Servicing date exceeded. Servicing to be scheduled.
	Hourglass steady	No need for servicing.

For AIR ONE PLUS® / AXTAIR AUTOMORPHO AUTO PLUS® /AXTAIR XXL® compressors:  
The "spanner" icon LED lights on the compressor to indicate that servicing is required.

Spanner icon	Meaning
 FIXED LED	Means product servicing is required (17500 hours of operation).

Servicing compressors at least every 2 years of use is recommended. To facilitate the management and scheduling of this servicing, counters give information about compressor use and the time since last servicing. This information is available using the Aircare maintenance software.

**⚠ WARNING:**  
Servicing can only be carried out authorized persons, **contact your dealer or the manufacturer.**

**WEEE (WASTE ELECTRIC AND ELECTRONIC EQUIPMENT)**

In compliance with applicable regulations covering the elimination of Waste electric and electronic equipment (WEEE), Winncare France is committed to contributing to the processing and recycling of electric and electronic equipment at the end of its service life.

The holder must never dispose of electric and electronic equipment at the end of its service life in the environment or with household waste with other unsorted waste, so that it can be the subject of selective collection in order to be recycled or reused in order to protect the environment, protect natural resources and protect health.

**EQUIPMENT CONCERNED**

The devices and instruments that feature the symbol  (Waste electric and electronic equipment) meaning that at the end of their service life they must not be mixed with other types of waste, must be the subject of separate processing.

**BEFORE 13 AUGUST 2005: FOR ELECTRIC AND ELECTRONIC EQUIPMENT PUT ON THE MARKET**

Holders must finance and organize the elimination of electric and electronic equipment at the end of its service life themselves. Unless the equipment is replaced by an equivalent device or a device that has the same function. In that case, Winncare France applies the same waste electric and electronic equipment (WEEE) elimination procedure as for equipment put on the market after 13 August 2005.

**AFTER 13 AUGUST 2005: FOR ELECTRIC AND ELECTRONIC EQUIPMENT PUT ON THE MARKET**

Winncare France complies with applicable European regulations. Winncare France undertakes to finance and organize the elimination of such waste.


The procedures for the organization and elimination of WEEE are available from:

**WINNCARE France****4 LE PAS DU CHÂTEAU****85670 SAINT PAUL MONT PENIT****Tél. : +33(0)2 51 98 55 64****Fax : +33(0)2 51 98 59 07****email: [contact@winncare.fr](mailto:contact@winncare.fr) or on our website: [www.winncare.fr](http://www.winncare.fr)**

If the holder does not apply this procedure, they must recycle electric and electronic equipment at the end of its service life using an adapted method and at their own expense. They must communicate the necessary information to Winncare France. If the products are re-sold, the holder must make sure the new owner is perfectly informed of their duties in terms of recycling electric and electronic equipment at the end of its service life.

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**A3**


## List of references: VAXT/POMPE/ONEP

Rep.	Before		After **	
	References	References	References	References
1	PB2050521			Designation Self-adhesive isophonic foam
2	PB2050163			Power cord & clips
3	WS/01/0073*			Air distributor assembly
3.1	PB2050118*			Motor 0.83RPM CW+ connections
4	WS/01/0074*			Pump assembly
4.1	PB2050127*			Silentbloc
4.2	PB2050135*			Air pump
4.3	PB2050126*			Pump support
5	WS/01/0083*	WS/01/0084*		Front shell assembly
6	WS/01/0076*	WS/01/0087*		Rear shell assembly
6.1	PB2050171			Fixing hook
6.2	PB2050140	PB2050343		Fixing bracket hook
6.3	PB2050136			air filter
6.4	PB2050139			Filter cover
11	references depend on the language of your country***			
12	PB2050146	PB2050530		Rear Sticker
13	PB2050210*			Female Connector
14	WS/01/0081*			Electronic PCB+ support x 2
14.1	PB2050138*			PCB support
15	PB2050153*			230 volt female connector
16	PB2050154*			Ceramic Fuse 5x20
//	PB2050337*			Single-dose grease
//	VAXT/KITREV01*			Overhaul kit (all components required for the standard AirCare overhaul)

\* The use of these references requires the completion of AirCare training.

\*\* Compatible references if the serial number is **higher** than

- 18458140 (VAXT/POMPE/AUTOP)
- 18457072 (VAXT/POMPE/ONEP)
- 18455241 (VAXT6/POMPE/AUTO)

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## List of references: VAXT/POMPE/AUTOP

Rep.	Before References	After References	Designation
1	PB2050521		Self-adhesive isophonic foam
2	PB2050163		Power cord & clips
3	WS/01/0073*		Air distributor assembly
3.1	PB2050118*		Motor 0.83RPM CW+ connections
4	WS/01/0074*		Pump assembly
4.1	PB2050127*		Silentbloc
4.2	PB2050135*		Air pump
4.3	PB2050126*		Pump support
5	WS/01/0077*	WS/01/0085*	Front shell assembly
6	WS/01/0076*	WS/01/0087*	Rear shell assembly
6.1	PB2050171		Fixing hook
6.2	PB2050140	PB2050343	Fixing bracket hook
6.3	PB2050136		air filter
6.4	PB2050139		Filter cover
11	references depend on the language of your country		
12	PB2050146	PB2050530	Rear Sticker
13	PB2050210*		Female Connector
14	WS/01/0081*		Electronic PCB+ support x 2
14.1	PB2050138*		PCB support
15	PB2050153*		230 volt female connector
16	PB2050154		Ceramic Fuse 5x20
//	PB2050337*		Single-dose grease
//	VAXT/KITREV01*		Overhaul kit (all components required for the standard AirCare overhaul)


\* The use of these references requires the completion of AirCare training.

\*\* Compatible references if the serial number is **higher** than :

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
## List of references: VAXT/PAUTOP

Rep.	References	Designation
1	PB2050521	Self-adhesive isophonic foam
2	PB2050163	Power cord & clips
3	W/S/01/0073*	Air distributor assembly
3.1	PB2050118*	Motor 0.83RPM C/w+ connections
4	W/S/01/0074*	Pump assembly
4.1	PB2050127*	Silentbloc
4.2	PB2050135*	Air pump
4.3	PB2050126*	Pump support
5	W/S/01/0085*	Front shell assembly
6	W/S/01/0087*	Rear shell assembly
6.1	PB2050171	Fixing hook
6.2	PB2050343	Fixing bracket hook
6.3	PB2050136	air filter
6.4	PB2050139	Filter cover
11	references depend on the	Ax/air side sticker
12	PB2050530	Rear Sticker
13	PB2050532*	Female Connector
14	W/S/01/0081*	Electronic PCB+ support x 2
14.1	PB2050138*	PCB support
15	PB2050153*	230 volt female connector
16	PB2050154	Ceramic Fuse 5x20
17	PB2050553	QR CODE Sticker
//	PB2050337*	Single-dose grease
//	VAXTKITREV01*	Overhaul kit (all components required for the standard AirCare overhaul)

\* The use of these references requires the completion of AirCare training.

\*\*\*  
 PB2050566 = GB PB2050575 = IT PB2050584 = PT  
 PB2050569 = DE PB2050578 = NL PB2050587 = DK  
 PB2050572 = ES PB2050581 = PL PB2050590 = SK

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## List of references: VAXT/PONEP


Rep.	References	Designation
1	PB2050521	Self-adhesive isophonic foam
2	PB2050163	Power cord & clips
3	WS/01/0073*	Air distributor assembly
3.1	PB2050118*	Motor 0.83RPM CW+ connections
4	WS/01/0074*	Pump assembly
4.1	PB2050127*	Silentbloc
4.2	PB2050135*	Air pump
4.3	PB2050126*	Pump support
5	WS/01/0084*	Front shell assembly
6	WS/01/0087*	Rear shell assembly
6.1	PB2050171	Fixing hook
6.2	PB2050343	Fixing bracket hook
6.3	PB2050136	air filter
6.4	PB2050139	Filter cover
11	references depend on the language of your country***	Side Sticker Axtair
12	PB20503530	Rear Sticker
13	PB2050532*	Female Connector
14	WS/01/0081*	Electronic PCB+ support x 2
14.1	PB2050138*	PCB support
15	PB2050153*	230 volt female connector
16	PB2050154*	Ceramic Fuse 5x20
16	PB2050553	Sticker QR CODE
//	PB2050337*	Single-dose grease
//	VAXT/KITREV01*	Overhaul kit (all components required for the standard AirCare overhaul)

\* The use of these references requires the completion of AirCare training.

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Rep.	Before References	After ** References	Designation
1	PB2050521		Self-adhesive isophonic foam
2	PB2050163		Power cord & clips
3	WS/01/0073*		Air distributor assembly
3.1	PB2050118*		Motor 0.83RPM CW+ connections
4	WS/01/0074*		Pump assembly
4.1	PB2050127*		Stentbloc
4.2	PB2050135*		Air pump
4.3	PB2050126*		Pump support
5	WS/01/0075*	WS/01/0086**	Front shell assembly
6	WS/01/0076*	WS/01/0087**	Rear shell assembly
6.1	PB2050171		Fixing hook
6.2	PB2050140*	PB2050343**	Fixing bracket hook
6.3	PB2050136		air filter
6.4	PB2050139		Filter cover
11			references depend on the language of your country***
12	PB2050146*	PB2050530**	Rear sticker
13	WS/01/0080*	WS/01/0080**	Female air/data connector + Axensor Bluetooth card
14	WS/01/0081*		Electronic PCB+ support x 2
14.1	PB2050138*		PCB support
15	PB2050153*		230 volt female connector
16	PB2050154		Ceramic Fuse 5x20
17	PB2050178*		Bluetooth PCB support
18	PB2050553		QR CODE Sticker
//	PB2050337*		Single-dose grease
//	VAXT/KITREV01*		Overhaul kit (all components required for the standard AirCare overhaul)

\* The use of these references requires the completion of AirCare training.

\*\* Compatible references if the serial number is **higher** than :

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
Rep.	References	Designation
1	PB2050521	Self-adhesive isophonic foam
2	PB2050163	Power cord & clips
3	WS/01/0073*	Air distributor assembly
3.1	PB2050118*	Motor 0.83RPM CW+ connections
4	WS/01/0074*	Pump assembly
4.1	PB2050127*	Silentbloc
4.2	PB2050135*	Air pump
4.3	PB2050126*	Pump support
5	WS/01/0092*	Front shell assembly
6	WS/01/0087*	Rear shell assembly
6.1	PB2050171	Fixing hook
6.2	PB2050343	Fixing bracket hook
6.3	PB2050136	air filter
6.4	PB2050139	Filter cover
11	references depend on the language of your country***	Axtair side sticker
12	PB2050530	Rear sticker
13	WS/01/0093*	Female air/data connector + wiring
14	WS/01/0081*	Electronic PCB+ support x 2
14.1	PB2050138*	PCB support
15	PB2050153*	230 volt female connector
16	PB2050154	Ceramic Fuse 5x20
17	PB2050553	QR CODE Sticker
//	PB2050337*	Single-dose grease
//	VAXT/KITREV01*	Overhaul kit (all components required for the standard AirCare overhaul)


\* The use of these references requires the completion of AirCare training.

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## List of references: VAXT/POMPE/ XXL

Rep.	References	Designation
1	PB2050521	Self-adhesive isophonic foam
2	PB2050163	Power cord & clips
3	WS/01/0073*	Air distributor assembly
3.1	PB2050118*	Motor 0.83RPM CW+ connections
4	WS/01/0074*	Pump assembly
4.1	PB2050127*	Silentbloc
4.2	PB2050135*	Air pump
4.3	PB2050126*	Pump support
5	WS/01/0088*	Front shell assembly
6	WS/01/0087*	Rear shell assembly
6.1	PB2050171	Fixing hook
6.2	PB2050343	Fixing bracket hook
6.3	PB2050136	Air filter
6.4	PB2050139	Filter cover
11	references depend on the language of your country***	Side sticker Axtair
12	PB2050530	Rear sticker
13	PB2050210*	Female Connector
14	WS/01/0081*	Electronic PCB+ support x 2
14.1	PB2050138*	PCB support
15	PB2050153*	230 volt female connector
16	PB2050154*	Ceramic Fuse 5x20
17	PB2050553	QR CODE Sticker
//	PB2050337*	Single-dose grease
//	VAXT/KITREV01*	Overhaul kit (all components required for the standard AirCare overhaul)






\* the use of these references requires the AirCare training

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